
ENG

Huvitz

CORNEAL TOPOGRAPHY **HTG-1**

USER MANUAL

IMPORTANT NOTICE

Electromagnetic waves caused by devices such as mobile phones, transceivers, and radio-controlled toys may lead to malfunctions of this product. Please keep objects that would affect the product away from it.

The information in this user manual has been carefully checked, and is believed to be entirely accurate at the time of publication. However, HUVITZ assumes no responsibility for any potential errors or omissions resulting from the use of the information contained in this manual.

HUVITZ reserves the right to make changes to this product or its specifications at any time without particular notice and is not required to update this documentation to reflect such changes.

■ Revision History

Revision	Date	Approval	Description
A	2021.05.27	 박 상 준	First issued
B	2023.04.21	 박 상 준	Added a NOTE for the Automatic Tracking function Added a NOTE regarding errors in measurement
C	2024.06.03	 채 기 병	Updated the information of U.S. representative address
D	2024.08.21	 박 상 준	Added the Dry Eye function Updated GUI main screens
E	2025.07.04	 박 상 준	Updated modified GUI main screens Updated modified headrest image and dimension Updated GUI screen shoot Changed exterior design (reduced body size) Changes to the final version of the Dry Eye test result screen screenshot (Comparison function added)

9000ENG0095-E
(2025/07/04)

© 2021 HUVITZ Co., Ltd.

38, Burim-ro 170beon-gil, Dongan-gu, Anyang-si, Gyeonggi-do, 14055, Republic of Korea

All rights reserved.

Under copyright laws, this manual may not be copied, in whole or in part, without the prior written consent of HUVITZ Co., Ltd.

CONTENTS

SAFETY PRECAUTIONS	4
1.1. Overview.....	4
Symbol Information.....	5
2.1. Usage Precautions.....	8
2.2. Environmental Considerations	10
2.3. Safety Warnings.....	12
2.4. Safety Cautions	17
INTRODUCTION.....	18
3.1. System Outline	18
3.2. Intended Use	18
3.3. Classifications	18
3.4. Contraindications	18
3.5. Patient Requirements	19
3.6. Operating Principles.....	19
3.7. Applied Standard List.....	19
SYSTEM OVERVIEW	20
4.1. Configuration and Functions	20
4.2. Main Body Screen Description	22
INSTALLATION PROCEDURE.....	39
5.1. System Installation.....	39
Operation	41
6.1. Software	41
6.2. SETUP Mode.....	43
6.3. DICOM SETUP Mode.....	47
6.4. Device Inspection	50
6.5. General Operation	51
6.6. MAINTENANCE.....	55
6.6.1. After operation.....	55
6.6.2. Cleaning and Instrument Protection.....	56
6.7. Send Data to Web-Viewer.....	57
TROUBLESHOOTING GUIDE	59
Specifications and Accessories.....	60
8.1. Standard Accessories	60
8.2. Specifications.....	61
8.3. Drawings of System.....	61
EMC INFORMATION	62
SERVICE INFORMATION.....	64
SOFTWARE LICENSE AGREEMENTS	65

1

SAFETY PRECAUTIONS

1.1. Overview

Safety is everyone's responsibility. The safe use of this instrument is largely dependent upon the installers, users, operators, and managers. It is prerequisite to read and understand these specifications before installing, using, cleaning, fixing or revising. Fully understanding the whole instructions must be the first priority. For this reason, the following safety notices have been placed appropriately within the text of this manual to highlight safety related information or information requiring special emphasis. All users, operators, and maintainers must be familiar with and pay particular attention to all signs of Warnings and Cautions.

WARNING

"Warning" indicates the presence of a hazard that could result in severe personal injury, death, or substantial property damage if ignored.

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CAUTION

"Caution" indicates the presence of a hazard that could result in minor injury or property damage if ignored.

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NOTE

This is used to emphasize essential information.

Be sure to read this information to avoid operating the device incorrectly.

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


2

Symbol Information

The International Electrotechnical Commission (IEC) has established a set of symbols for medical electronic equipment which classify a connection or warn of any potential hazards. The classifications and symbols are shown below.

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human health. For more detailed information about disposal of your old appliance, please contact your city office, waste disposal service or the shop where you purchased the product.
(Symbole WEEE- EU seulement)

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Refer to instruction manual/booklet
(Se reporter au manuel d'instructions / brochure)

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MEDICAL – 355544 EQUIPMENT
AS TO ELECTRICAL SHOCK, FIRE AND MECHANICAL HAZARDS ONLY IN
ACCORDANCE WITH ANSI/AAMI ES60601-1:2005/A2:2021, CAN/CSA-C22.2 No. 60601-1
(Amendment 2:2022)

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2.1. Usage Precautions

This instrument has been developed and tested in conformity with domestic and international safety standards and regulations. This ensures that the instrument maintains a high degree of safety. By law, the manufacturer must thoroughly inform the user about the safety aspects of the instrument. It is imperative that the user should use this instrument according to the manual for their own safety. Therefore, please read the instructions in the manual carefully and fully understand them before turning on the instrument. For more detailed information, please contact our Customer Service Department or one of our authorized representatives.

CAUTION

For use of equipment in rated voltage less than 125Vac, minimum 6A, Type SJT or SVT, 18/3AWG, 10A, max 3.0m long: One end with Hospital Grade Type, NEMA 5-15P Other end with appliance coupler.

For use of equipment in rated voltage less than 250Vac, minimum 6A, Type SJT or SVT, 18/3AWG, 10A, max 3.0m long: One end terminated with blade attachment plug(HAR) Type, NEMA 6-15P.

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
WARNING

Use instrument that comply with IEC60601-1 in the patient environment. [The figure below show]

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If an instrument that does not comply with IEC 60601-1 is to be used, use an isolation transformer.

If a person handling a conductive part of the system comes into contact with a patient at the same time, hazard may occur due to leakage current exceeding the value specified in the applicable standard. Be careful not to touch patients when connecting or removing the power plug or cable connectors.

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 CAUTION

This instrument includes lithium battery. This hazardous material needs to be disposed of properly to limit environmental pollution. Please contact to the professional waste disposal company.

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 CAUTION

Do not install any software on device without our consent.

The manufacturer is not responsible for any failure due to random installation.

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 CAUTION

1 Light Hazard Group I (according to the ISO15004-2)

This device is classified as a Group I. The LED used for the device is safe under expected use.

However, observe the following precautions when using the device.

- 1.1 Do not direct LED beams to human eyes when unnecessary.
- 1.2 Do not look into the objective lens for a prolonged time.
- 1.3 If the device cannot be used appropriately or any problem cannot be solved, Stop using the device immediately, and contact Huvitz or your authorized distributor.

 NOTE

The device must be operated only by a, or under direct supervision of properly trained and qualified person/s such as ophthalmologist.

a. Age: adult (> 21 years of age)

b. Occupation: Doctor, ophthalmologist

c. Level of training: Recommend completing training how to use it is.

d. Education: Basic training is required for the use of equipment.


(Obtains a certificate of qualification such as optics, medical school, etc., after completing a certain amount of education.)

2.2. Environmental Considerations


Avoid the following environments for operation or storage:

Places where the instrument is directly exposed to moisture.
Do not operate the instrument with wet hands.

Places where the instrument is directly exposed to sunlight


Places where the instrument is exposed to ultraviolet rays

Places where there are large changes in the temperature


The optimal temperature range for normal operation is 10°C to 35°C, and the humidity range is 30% to 90%.

Places where appliances that generate heat are in close proximity

Places with high humidity and heat emission issues

Places where the instrument is exposed to excessive shock or vibration

Places where the instrument is exposed to chemical materials or explosive gas

Be careful not to let any dust or metal get inside the instrument.

Do not disassemble or open the instrument. The manufacturer does not take responsibility for any problems caused by such actions.

Do not block the heat vent on top of the instrument.

Do not plug the AC power cable unless all instrument parts are completely connected. Otherwise, the instrument will be damaged.

Do not pull out the power cable by the cord.

This instrument can withstand the following conditions:

1. Operation

- An ambient temperature range of 10°C ~ 35°C (50°F ~ 95°F)
- A relative humidity range of 30% ~ 90% (with non-condensing)
- An atmospheric pressure range of 800 ~ 1060hpa
- An illuminance of the measuring room at 10~100 lux

2. Transportation

- An ambient temperature range of -40°C ~ 70°C (-40°F ~ 158°F)
- A relative humidity range of 10% ~ 95%
- An atmosphere pressure range of 500 ~ 1060hpa

3. Storage

- An ambient temperature range of -10°C ~ 55°C (14°F ~ 131°F)
- A relative humidity range of 10% ~ 95% (with non-condensing)
- An atmosphere pressure range of 700 ~ 1060hpa

Avoid environments where the device is exposed to excessive shocks or vibrations.

Do not expose products or packing to environmental conditions outside of those specified above.


2.3. Safety Warnings

WARNING

1. This is an electric medical device. Use is limited to doctors or persons qualified by the law of each country.
2. Do not make a diagnosis based on a single captured image. Doctors are responsible for making the final diagnosis based on the present and past medical records of the patient such as captured images. Without sufficient information, proper diagnosis may not be made.
3. This device must not be used in an area that is in danger of explosions and in the presence of flammable, explosive, or volatile solvent such as alcohol, benzene or similar chemicals.
4. Do not place or store this instrument in humid area. Do not expose the device to water splashes, dripping water, or sprayed water. Do not place containers with fluids, liquids, or gases on top of this instrument.
5. The device must be operated by a trained and qualified person or under his or her supervision.
6. Repair of this instrument must be conducted by HUVITZ's service technicians or other authorized persons.
7. When a user takes care of its maintenance, the user must follow the instructions in the user and service manuals. Any additional maintenance work must be conducted by a HUVITZ service technician or an equivalently qualified person.
8. Manufacturers are not responsible for the damages caused by unauthorized alterations. Such tampering will forfeit any rights to receive services during the term of guarantee.
9. This instrument should be used with accessories supplied by HUVITZ. If the consumer is to use accessories from another manufacturer, their safety and usability must be checked and proved by HUVITZ or the accessories' manufacturer.
10. Only those who have undergone proper training or education may install, operate, and maintain this instrument.
11. Do not apply excessive force to cable connections. If you find it difficult to connect the cable, check if the connector (plug) is appropriate for the socket. If you have damaged the connector or the socket, you must contact a qualified service technician for repair..
12. Do not pull the instrument's cables. Hold the plug and pull it out when disconnecting a cable.
13. Do not block any ventilation outlet necessary for proper heat dissipation.
14. If smoke, sparks or any abnormal noise or smell is noticed coming from the instrument, please switch the power off immediately and pull out the plug.
15. To avoid the risk of electric shock, this instrument must only be connected to power supply with a protective ground.

16. Do not When placing the instrument, do not place it in a location that makes it difficult to connect or disconnect cables (power cable, external device-connecting cable, etc.).
17. If you intend to connect an external device to an input/output signal or other connectors, they must meet the relevant IEC standards (IEC 60950 for IT devices and IEC 60601-1 series for medical electrical devices). In addition, all such combination-system-shall comply with the standard IEC60601-1 harmonized national standard or the combination. If, in doubt, contact qualified technician or your local representative. The operator should not touch the patient and accessible male parts of the SIP/SOP connectors simultaneously.
18. When carrying the instrument, hold its bottom left (A) and right (B) sides.

19. Do not touch the instrument directly if the user has a wound in their hand or is allergic to the material used for the contact area.


Part Name	Material
LCD Touch pannel	Glass
Joystick / button	ABS + Silicon, Aluminum(A6061 T6) + Anodizing
Power switch	PC + PA66
Cover	ABS
Chin Rest	ABS

20. Refrain from using the instrument on a patient who is sensitive to light. (e.g. photophobia)
21. When instrument is send back to A/S center for repair or maintenance, or before authorized service man is arrived at the place for repair or maintenance, wipe the surfaces of the instrument (especially, the parts that come into contact with the patient) with a clean cloth dampened with rubbing alcohol.
22. In the event of a serious incident involving the device, the user shall report it to the manufacturer and the competent authority of the Member State in which the user and/or patient are established.
23. If a user uses power save supported by Windows 10 except for the power save in User setting, it causes some trouble in HTG-1. The manufacturer isn't responsible for the problem.
24. User must not change the setting supported by the manufacturer. This change might make some trouble in HTG-1. The manufacturer isn't responsible for the problem.
25. Immediately after the test, you may not be able to see clearly due to glare. It is recommended to sit back, relax and move around.
26. When packing box is unintentionally open before use or damaged, please call A/S center.
27. Make sure that you always have an up-to-date virus scanner and firewall installed on your system.
28. Install and operate software within secure operating environment that allows only authorized users to access and that the system network is equipped with Window firewall built-in Windows system, windows Defender antispysware tools and other common security tools /application systems
29. Avoid resetting or turning off the PC when the HTG-1 running.

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INTRODUCTION

3.1. System Outline

The Huvitz Corneal Topographer HTG-1 is a device for measuring and analyzing the cornea topography. The Corneal Topographer HTG-1 optically measures eye components such as Corneal curvature radius, Corneal cylinder axis, White to White, Pupil diameter and dry eye screening.

HTG-1 software based on Windows 10 operating system provides various tests of the anterior corneal surface. Along with mapping of the corneal surface, HTG-1 provides various analyses of the patient's corneal surface.

HTG-1 Software Specifications:

- Database management
- Corneal wavefront analysis (Zernike)
- Contact lens fitting software
- On-board fluorescein acquisition
- Auto-focus system
- Corneal analysis
- Integrated pupillometer
- Keratoconus screening
- Report management
- Reviewing previous patient examinations
- Dry eye screening

3.2. Intended Use

The Corneal Topographer HTG-1 is corneal analyzer as acquiring images of the cornea and measuring corneal topography. The rings of the placido reflected by the illuminated cone are used to geometrically calculate the topographic map of the cornea in the image.

3.3. Classifications

- Classification of product: Class IIa according to Medical Device Regulation 2017/745(MDR) Annex VIII Rule 10
- Resistance against electric shock: Class I (earthed)
- Protection class against electric: Type B (Headrest, chinrest paper)
- Classification of Light hazard: Group I (EN/ISO15004-2 Standard)

3.4. Contraindications

This device should not be used for:

- Patients who are hypersensitive to light.
- Patients who recently underwent photodynamic therapy
- Patients taking medication that causes photosensitivity

There is no any known residual risk other than risks arising from violating contraindication. But, when unexpected event occurs or you encounter different problems, please contact your local distributor or Authorized Representatives.

3.5. Patient Requirements

Patients who are examined with this instrument must maintain concentration for a few minutes and follow the instructions below.

- Firmly position your face on the chin rest and headrest.
- Keep your eyes open.
- Understand and follow the given instructions during the examination.

If patients do not adhere to the instructions above, proper images cannot be obtained.

3.6. Operating Principles

Corneal topography maps the corneal surface, a key refractive element of the eye. It uses a Placido disc pattern, reflected via a 50% beam splitter attached to the body of the microscope. Advanced image processing algorithms then perform edge detection on the reflected pattern to analyze corneal curvature.

3.7. Applied Standard List


- IEC/EN 60601-1: MEDICAL ELECTRICAL EQUIPMENT
 - Part 1: General requirements for safety
- IEC/EN 60601-1-2: Medical electrical equipment Part 1: General requirements for safety
 - Collateral Standard: Electromagnetic compatibility-requirements and tests
- ISO 15004-1: Ophthalmic instruments
 - Fundamental requirements and test methods
 - General requirements applicable to all ophthalmic instrument
- ISO 15004-2: Ophthalmic instruments - Fundamental requirements and test methods
 - Part 2: Light hazard protection
- ISO 19980: Ophthalmic instruments – Corneal topographers

4

SYSTEM OVERVIEW


4.1. Configuration and Functions

■ Front

No.	Part	Name	Description
1	Display	LCD	A monitor that displays the image to be captured and the icons that control the instrument, it can be tilted and adjusted according to the user's eye level.
2	Body	Chin rest button	A button that moves the chin rest up and down.
3		Joystick	A tool used for moving the instrument back-forth/left-right/up-down to align it with the patient's eye. After finishing the alignment, the user can begin measurement.
4		User lock	Fixes the instrument's body to its base.
5		Heat vent	A window where the internal heat is released outside.
6	Base	Power button	A button for powering the instrument's PC ON/OFF. When the power is ON, the button is lit white.
7		Power switch	A button for powering the instrument ON/OFF.
8		Power inlet	An inlet for connecting a power cable to the instrument.
9	Headrest	Eye level mark	Mark for indicating base height of patient's eye.

Rear


No.	Part	Name	Description
1	Headrest	Headrest	Part where the patient places their forehead
2		Chin rest	Part where the patient places their chin
3	Body	Placido's disk	Enables to map the curvature of the cornea.
4		Focus LED	LEDs to examine patient's eye.
5	Base	USB port	A set of ports for USB communication (3 ports)
6		LAN port	A port for LAN communication (1 port)
7		HDMI port	A port for connecting an external display device (1 port)
8		Foot switch port	A port for connecting a foot switch (1 port)

NOTE

Interference between USB 3.0 and 2.4 GHz wireless signals may cause USB connectivity and data transfer speed issues. To prevent such issues, use an extension cable or change the wireless channel on the Wi-Fi device to reduce the interference. In addition, using only USB 2.0 devices is another possible solution since this can also reduce interference. If the problem persists, consider using a wired connection or upgrading the connection to a 5 GHz wireless network.


Bottom

No.	Part	Name	Description
1	Base	Packing lock	Fixes the instrument's body and base during transportation.

4.2. Main Body Screen Description


MEASURE Screen

No.	Name	Function
1	Patient Information	Displays patient ID and name in examination.
2	Patient	Goes back to the patient list screen.
3	Result	Goes to the result screen. The result screen shows detailed information on the measured data.
4	AUTO S	Shows auto shooting selection status.
5	AUTO T	Shows auto tracking selection status.
6	Setup	Changes the measurement settings.
7	OD&OS Info.	Displays information on the position (left/right) of the eye that is currently being measured.
8	Camera Image	Displays the real-time camera image.
9	Measurement Value	Displays the most recent measurement results. - K1: Flat Keratometry - K2: Steep Keratometry - Axis: Angle degree of K1 - Cyl: Cylinder value - WTW: White to white
10	Measurement Mode	Changes the measurement mode.
11	Measurement Scale	Changes the measurement unit(Diopter <-> Radius)

RESULT Screen

1. Keratometry result screen (KER)


No.	Name	Function
1	Map/Profile	Map: Displays the topography map with various options. Profile: Displays the curvature along the meridian.
2	Topo map image	Displays the topo-map image according to the set options.
3	Eye/Map/Ring	Enables the option to turn ON/OFF the display of eye, topo-map, and topo-ring images
4	Axial/Tangential/ Refractive/Elevation	Provides options for a topo-map for selection.
5	Zoom	Displays the map in full-screen.
6	Topography Result	Displays the detailed values (SimK, Meridian, Keratoconus) based on the selection. SimK: Displays the measured SimK data. Meridian: Displays the measured data for the area ranged 3 to 7 mm or 2 to 6 mm from the center point. (K1: Flat Keratometry, K2: Steep Keratometry, Avg: Average of Keratometry value in the zones, Cyl: Cylinder and Angle degree of K1) Keratoconus: Displays the expected results for keratoconus. (Non-Keratoconus, Keratoconus Suspicious, Central or Peripheral Steepening Keratoconus).
7	Scale	Set the topo-map options, such as color scale, step, and overlay.
8	Display Type	Provides selection options for Full Image, 4 Images, and Comparison. Full Image: Displays one map in an enlarged version. 4 Images: Displays all maps for Axial, Tangential, Refractive Power, and Elevation options on one screen. Comparison: Compares the previously measured map data with the current data and displays the comparison map together.
9	Eye Image & Meridian	Displays an eye image marked with steep and flat meridians.
10	Control Meridian	Changing a meridian's angle updates the corresponding graph.
11	Value Graph & Differential Graph	Shows the radius/diopter based on the meridian types. On the differential graph, you can check the difference between the steepest meridian and the flattest meridian.

2. Zernike result screen (ZER)


No.	Name	Function
1	Map/PSF	Map: Displays the aberration map. PSF: Displays the point spread image and the applied chart image.
2	Low Order Aberrations/ High Order Aberrations	Displays the maps of low-order/high-order aberrations.
3	Expansion Coefficients	Displays the Zernike coefficient values in a graph.
4	Point Spread Function	Displays a PSF image consisting of Zernike coefficients.
5	Chart Images	Displays an image where the PSF is applied to the selected chart.

3. Pupillometry result screen (PUP)


No.	Name	Function
1	Pupil/Graphs	Pupil: Displays continuous images of the pupil. Pupil changes can be observed. Graph: Displays the latency and decentration graph.
2	Captured Image	Displays the image captured during measurement.
3	Display Options	Among the options below, displays the selections together on the screen.
4	Measurement Values	Displays the measured data.
5	Latency	Displays the measured pupil size in a latency graph.
6	Decentration	Displays how far the eye's center and the pupil's center are out of alignment.

4. White to White result screen (WTW)

No.	Name	Function
1	WTW Bar	Displays the result of the WTW measurement.
2	Edit	Edits the WTW bar status on the edit window.
3	Values	Displays the White to White value.
4	Edit Window - Edit Bar	The bar turns pink when clicked, which indicates that the edit mode is now active. While the bar is active, you can move it to change the WTW position.
5	Edit Window - Detailed Move Buttons	Since it is difficult to make a detailed move for the bar on a small screen, these buttons are used for fine adjustments.
6	Edit Window - Reset	Relocates the bar to the original position.
7	Edit Window - Save and Close	After finishing the editing, saves and closes the edit window.

5. Contact Lens Fitting result screen

No.	Name	Function
1	Lens Fitting Map/ Corneal Curvature Map	Displays the lens fitting map or corneal curvature map of the subject's selected lens.
2	Select Corneal Curvature Map	Provides options for a corneal curvature map for selection. (Axial / Tangential Maps are available)
3	Select Lens Fitting Map	Displays the lens fitting map if there is a selected lens.
4	Major/Minor Axis Lens Fitting Graph	Displays the distance between the cornea and the lens for the major/minor axes for the selected lens.
5	Select Lens Manufacturer	Provides options for a lens manufacturer for selection.
6	Select Lens	After selecting a manufacturer, selects a lens from the list of lenses manufactured by them.
7	Modify Lens Parameters	Modifies the parameters of the selected lens.
8	Corneal Curvature Information	Displays the information on the subject's corneal curvature.

5.1. Lens Parameter Modification Screen

Lens Parameters

Base Curve

Curvature Radius(r0) ①
- 8.25 +

Diameter(mm) ②
- 9.8 +

Curvature Diff(dr0) ③
- 0 +

Eccentricity ④
- 0.9 +

Peripheral Curve 1 ⑤

Curvature Radius(r1)
- 0 +

Width(mm)
- 0 +

Curvature Diff(dr1)
- 0 +

Eccentricity
- 0 +

Peripheral Curve 2 ⑥

Curvature Radius(r2)
- 0 +

Width(mm)
- 0 +

Curvature Diff(dr2)
- 0 +

Eccentricity
- 0 +

Peripheral Curve 3 ⑦

Curvature Radius(r3)
- 0 +

Width(mm)
- 0 +

Curvature Diff(dr3)
- 0 +

Eccentricity
- 0 +

Default ⑩
Cancel ⑧
Apply ⑨

No.	Name	Function
1	Set Lens Curvature	Set the curvature of the lens.
2	Set Lens Length (diameter/width)	Set the diameter/width of the lens. Base Curve Diameter (mm): Indicates the full length of the lenses. Peripheral Curve Width (mm): Indicates the length of the corresponding peripheral curve.
3	Set Curvature Difference	Set the difference in the curvature of the lens.
4	Set Eccentricity	Set the eccentricity of the lens.
5	Set Primary Curve	Set the curvature, length and the difference/eccentricity of the curvature for the primary curve.
6	Set Secondary Curve	Set the curvature, length and the difference/eccentricity of the curvature for the secondary curve.
7	Set Tertiary Curve	Set the curvature, length and the difference/eccentricity of the curvature for the tertiary curve.
8	Cancel Settings	Cancels the current settings and returns to the lens fitting screen.
9	Apply Settings	Saves the current settings and returns to the lens fitting screen.
10	Default Setting	Automatically sets the lens setting that most closely match the currently measured result.

6. Dry Eye Screening result screen
6.1. Questionnaire(OSDI)

QUESTIONNAIRE (OSDI)

Have you EXPERIENCED any of the following during the last week: 1

	All of the time	Most of the time	Half of the time	Some of the time	None of the time
1. Eyes that are sensitive to light?	4	3	2	1	0
2. Eyes that feel gritty?	4	3	2	1	0
3. Painful or sore eyes?	4	3	2	1	0
4. Blurred vision?	4	3	2	1	0
5. Poor vision?	4	3	2	1	0

< 1 / 3 >

QUESTIONNAIRE (OSDI)

Have PROBLEMS WITH your eyes limited you in performing any of the following during the last week:

	All of the time	Most of the time	Half of the time	Some of the time	None of the time
6. Reading?	4	3	2	1	0
7. Driving At Night?	4	3	2	1	0
8. Working With A Computer Or Bank Machine(ATM)?	4	3	2	1	0
9. Watching TV?	4	3	2	1	0

< 2 / 3 >

QUESTIONNAIRE (OSDI)


Have your eyes felt UNCOMFORTABLE in any of the following situations during the last week:

	All of the time	Most of the time	Half of the time	Some of the time	None of the time
10. Windy Conditions?	4	3	2	1	0
11. Places Or Areas With Low Humidity	4	3	2	1	0
12. Areas That Are Air Conditioned?	4	3	2	1	0

< 3 / 3 >


No.	Name	Function
1	Questionnaire Content	Fills out the questionnaire.

6.2. Kearatoscopy Image(Fluorescein)


No.	Name	Function
1	Type	Provides options to select the type of dry eye screening result.
2	Select capture image	Provides options to select the captured image to check the results.
3	Zoom	Displays the result on a bigger screen.

6.3. Lipid Layer screening


No.	Name	Function
1	Type	Provides options to select a dry eye screening result.
2	Select video recorded	Provides options to select the recorded video to check the results.
3	Play	Plays a video.
4	Zoom	Displays the result on a bigger screen and evaluate it.

6.4. Tear Meniscus Height screening


No.	Name	Function
1	Type	Provides options to select a dry eye screening result.
2	Select captured image	Provides options to select the captured image to conduct analysis.
3	Analysis Result	Displays the analysis result.
4	Detail Analysis	Displays the analysis to check and edit the selected image in detail.
5	Zoom	Zooms in/out the image for analysis.
6	Edit	Set and edit the selected area for analysis.
7	Comparison	Compares result with previous measurement data.

6.5. NIKBUT (tear film break-up time screening)


No.	Name	Function
1	Type	Provides options to select the desired dry eye screening result.
2	Select recorded video	Provides options to select the recorded video to check the results.
3	Analysis Result	Displays the analysis result.
4	Detail Analysis	Displays the analysis to check the selected video data in detail.
5	Color Map	Displays the broken area of the tear film in a map with the degree of the break-up.
6	Play	Plays a video or navigates a video frame by frame.
7	Eye & Map	Displays the eye/map images together on the video based on your selection.
8	Breakup Graph	Shows the break-up over time in a graph.
9	Comparison	Compares result with previous measurement data.

6.6. Meibography (meibomian gland screening)


No.	Name	Function
1	Type	Provides options to select a dry eye screening result.
2	Select captured image	Provides options to select a shot image to conduct meibography.
3	Analysis Result	Displays the analysis result.
4	Detail Analysis	Displays the analysis to check and edit the selected image in detail.
5	Area Mode	Provides “Grade” and “Edit” mode to evaluate meibography.
6	Detection Control	Set the detection level for meibography.
7	Edit	Set and edit the selected area for meibography.
8	Enhancement Enabled	Displays the enhanced image while highlighting.
9	Comparison	Compares result with previous measurement data.

6.7. Redness (redness degree screening)

No.	Name	Function
1	Type	Provides options to select a dry eye screening result.
2	Select image captured	Provides options to select a captured image to check the results.
3	Analysis Result	Displays the analysis result.
4	Zoom	Displays the result on a bigger screen.

6.8. Blink (blinking degree screening)


No.	Name	Function
1	Type	Provides options to select a dry eye screening result.
2	video List	Displays the list of recorded videos.
3	Analysis Result	Displays the analysis result.
4	Blink Graph	Displays the blinking moments over time in a graph.
5	Play	Plays the video.
6	Zoom	Displays the result on a bigger screen.

5


INSTALLATION PROCEDURE

5.1. System Installation

1. Place the main body on a stable table.
2. Loosen the two packing lock screw under the main body.

3. Unscrew the user lock lever on the body.

4. Attach the chin rest paper to the chin rest.
5. Connect external devices, if necessary.
 - (1) The PORT CAP can be opened by pushing the end of the bottom surface by hand.
 - (2) Connect mouse and keyboard, if necessary.
 - (3) Connect communication cable of external device.
 - (4) Close 'PORT CAP' with screw driver.
 - (5) The PORT CAP can be separated from the product and reassemble when connecting the cable. (Figure "A" below)




Figure "A"


6. Check if the power switch on the bottom right of the instrument's base is turned OFF (the O part is pressed).
7. Connect the power cable to the power inlet, then connect the other side to an AC outlet.
(The power cable is included in 8.1. Accessories)

8. If any external devices are connected, turn them on first.
9. Press the I part of the power switch on the bottom right of the instrument's base to supply power.
10. Press the power button on the bottom front of the instrument's base to power the instrument's PC.

11. Check if any error occurs during the initialization process.
Wait and do not touch anything until the initialization is completed.
12. Use the joystick and chin rest button on the body to check if there is any problem with the body's back-forth/left-right/up-down movements and the chin rest's up-down movement.


6

Operation

6.1. Software

1. Double-click the program icon on the Desktop.
2. Once the initialization process is finished, enter the user ID and password to log in.
3. Press the patient icon () and fill in the patient's information. If the patient is already registered, skip this step.

New Patient

Patient ID:

First Name:


Last Name:

Birth: 19840730

Gender:

Physician:

4. Check if the entered patient information is correct, then select a newly registered patient or an already registered patient from the list.

PATIENT LIST

Search:

Icons:     

ID	Name	Gender	Age	Update
☒ 00002	Pupillometry	Male	1983-08-17	2024-04-09
☒ 00001	Myopia Management Sample	Male	2006-07-31	2024-03-14
☒ 00006	Patient DCM1	Male	1983-08-17	2024-03-07
☒ 00004	Keratoconus Sample	Male	1983-08-17	2024-02-19

Left sidebar icons:     

5. To send the patient information to a web viewer or delete it, select the left circle of a desired patient ID and click the TRANSFER () button or the DELETE () button.

6. Once you select a patient, the screen changes as shown below.

7. Click the Measure icon () to enter the measurement mode.

6.2. SETUP Mode


[Select the section to setup]

- The setup consists of six sections: System, Patient, Measure, Connectivity, Report, Information. Select the desired section to setup displayed on the left.
- You can navigate through the different pages within the sections using the arrow icons [< , >]

[How to change the settings]

- Select the item you wish to change from the screen.
- If you press the OK button after changing the settings automatically, it saves any changes and exits the SETUP screen.
- If you press the CANCEL button, it exits the SETUP screen without saving the changes.

1. System Settings

Language	Set the language.
Device Name	Set the device name.
Sleep Time	Set Sleep mode time.
Touch Keyboard	Set whether to turn ON/OFF the touch keyboard function for use.
Login Page	Set whether to turn the login page ON/OFF for use.
Contact Lens	Imports/exports contact lens data.

2. Patient Settings

SETUP

System

Patient

Measure

Connectivity

Report

Information

Patient List Size

50

100

150

200

500

Today List

On

Off

PID Prefix

Reset

PID Postfix

Reset

PID Number Length

5

6

7

8

Data Format

YMD

MDY

DMY

Staff Management

Physicians

Default

<

1 / 1

>

Cancel

OK

Patient List Size	Set the number of patients displayed on the patient list.
Today List	Set whether to turn the Today List (list of patients who visited today) ON/OFF for use.
PID Prefix	Set the prefix of a patient ID.
PID Postfix	Set the postfix of a patient ID.
PID Number Length	Set the length of a patient ID.
Date Format	Set the date format (year, month, day).
Staff Management	Register staff.

3. Measure Settings

SETUP

System

Patient

Measure

Connectivity

Report

Information

Auto Tracking

On

Off

Auto Shoot

On

Off

Map Type

Axial

Tangential

REF Power

Elevation

Topo Recovery

On

Off

Measurement Scale

Diopter

Millimeter

Keratometry Display Area

2-4-6mm

3-5-7mm

Cylinder Notation

Negative

Positive

Default


<

1 / 2

>


Cancel

OK


Auto Tracking	Set whether to turn ON/OFF the auto tracking function.
Auto Shoot	Set whether to turn ON/OFF the auto measurement function.
Map Type	Select which map type to use as default (Axial / Tangential / Refractive Power / Elevation).
Topo Recovery	Set whether to turn ON/OFF the topo-map recovery function for use.
Measurement Scale	Select the default unit for measuring (Diopter / Millimeters).
Keratometry Display Area	3-5-7 mm / 2-4-6 mm
Cylinder Notation	Positive / Negative
Asphericity	e / SF / p / Q

4. Connectivity Settings

Auto Data Transfer	Set whether to turn ON/OFF the function that automatically sends the measured data to the web viewer upon saving
Web Viewer Server IP	Set the web viewer server IP address.
Web Viewer Sever Port	Set the web viewer server port.

5. Report Settings

Report Logo	Set the report logo image.
Organization	Assign the name of the organization that uses the system.
Auto Export	Set whether to turn ON/OFF the function that automatically creates a report file for the measured data upon saving for use.
Auto Export Folder	Set the path to automatically create a report file for the measured data upon saving.
Export Topo Map	Set the topo-map to be used for the report.

6.3. DICOM SETUP Mode

1. Worklist Settings

SETUP

System

Patient

Measure

Connectivity

Report

Information

Dicom

Worklist

Server IP

Server Port

Server AE

Client AE

Reset

Reset

Reset

Reset

Default

<

1 / 5

>

Cancel

OK

Server IP	Set the IP address of the PC where the server program is installed.
Server Port	Set the server port number.
Server AE	Set the server AE.
Client AE	Set the client AE.

2. MPPS Settings

SETUP

System

Patient

Measure

Connectivity

Report

Information

Dicom

MPPS

Server IP

Server Port

Server AE

Client AE

Reset

Reset

Reset

Reset

Default

<

2 / 5

>

Cancel

OK

Server IP	Set the IP address of the PC where the server program is installed.
Server Port	Set the server port number.
Server AE	Set the server AE.
Client AE	Set the client AE.

3. Storage Settings

SETUP

System

Patient

Measure

Connectivity

Report

Information

Dicom

Storage

Server IP

Server Port

Server AE

Client AE

Reset

Reset

Reset

Reset

Default

< 3 / 5 >

Cancel

OK

Server IP	Set the IP address of the PC where the server program is installed.
Server Port	Set the server port number.
Server AE	Set the server AE.
Client AE	Set the client AE.

4. General Settings

SETUP

System

Patient

Measure

Connectivity

Report

Information

Dicom

General

Manufacturer

Model Name

Institution

Station

Serial Number

Reset

Reset

Reset

Reset

Reset

Default

< 4 / 5 >

Cancel

OK

SETUP

System

Patient

Measure

Connectivity

Report

Information

Dicom

General

Set to Private UID

Private UID

Send MPPS

Transfer Syntax ▾

Modality ▾


5 / 5

Manufacturer	Set the name of the manufacturer.
Model Name	Set the model name.
Institution	Set the name of the institution.
Station	Set the name of the station.
Serial Number	Set the product's serial number.
Transfer Syntax	Set the transfer syntax.
Set to Private UID	Set whether to turn ON/OFF the Private UID.
Private UID	Set the private UID.
Send MPPS	Set whether to turn ON/OFF the MPPS send function.
Modality	Set the modality.


6.4. Device Inspection

1. Measurement using Model Eye


- (1) Remove the fixing pins from the chin rest, then remove the chin rest paper.

- (2) Mount the model eye illustrated below and fix it using the two fixing pins.

- (3) Using the chin rest button and joystick to align the center of the model eye lens with the center of the Placido ring.


- (4) Once you finish adjusting the focus, press the joystick button.


6.5. General Operation

1. Wipe the headrest and chin rest using a clean cotton swab or gauze, which is where a patient's forehead and chin comes into direct contact with the instrument.
When using the chin rest paper, remove one sheet from the top after every patient.

2. Use the joystick to align the left-right index mark (A) and the front-back index mark (B) on the body and base.


3. Instruct the patient to sit in front of the instrument.

4. Set the measurement mode.

Select one of the options displayed at the bottom of the screen: Kerato, Pupil, Dry Eye or Keratoscopy.

KERATO	Indicates that the Keratometer mode is ON.
PUPIL	Indicates that the Pupillometry mode is ON.
DRY EYE	Indicates that the Dry Eye mode is ON.
KERATOSCOPY	Indicates that the Keratoscopy mode is ON.

- In Keratometry mode, the Keratometry image is taken by pressing the button on the joystick, and the topography map and kerato value is calculated.
- In Pupillometry mode, the pupil image is taken by pressing the button on the joystick, and the pupil image and data is calculated.
- In Dry Eye mode, according to the selection of the desired mode (Lipid Layer, Tear Meniscus Height, NIBUT, Meibography, Blink, Redness).
- In Keratoscopy mode, either an image or video can be recorded for analysis and also supports Fluorescein screening

5. Select and set various functions from the left side of the screen.

(1) Set the Auto Shooting function.

	Indicates that the Auto Shooting function is ON.
	Indicates that the Auto Shooting function is OFF.

- If the Auto Shooting function is turned ON, the instrument automatically optimizes and takes the picture once the alignment and focus adjustment are completed.


(2) Set the Auto Tracking function.

	Indicates that the Auto Tracking function is ON.
	Indicates that the Auto Tracking function is OFF.


- If the Auto Tracking function is turned ON, the instrument automatically tracks the eye's center to align and adjust the focus when the patient's eye comes into the designated area.

6. Align the patient's eye to the eye-level mark on the headrest of the instrument.


- (1) Instruct the patient to put their chin on the chin rest.
- (2) Instruct the patient to press their forehead close against the headrest.
- (3) Adjust the height using the chin rest button on the body to align the patient's eye level with the eye-level mark on the headrest.

7. Instruct the patient to open their eyes wide and fix their gaze on the LED fixed inside.
8. Hold the joystick and adjust the body so that the patient's eye appears on the LCD screen.
9. Align and adjust the focus.
 - (1) Move the body up/down/left/right using the joystick until the center of the Placido ring (B) is aligned with the center of the cross (A). Once the two centers are aligned, check the focus indicator bar (C). Move the body back and forth using the joystick until the focus indicator bar (C) only shows yellow lines.


- (2) Adjust the joystick up-down/left-right/back-forth to align it until the target mark (yellow concentric circle, A) appears.

- (3) If the Auto Tracking function is turned ON, the instrument automatically aligns and adjusts the focus within the designated range.
- (4) If the body direction indicator (A) appears on the LCD screen during Auto Tracking, it means the moving

range has reached its limit. In this case, hold the joystick and move the body towards the indicated direction.

10. When in focus, press the joystick button to capture the image.

NOTE

The patient must keep their eyes wide open and not blink during the examination. Otherwise, the instrument cannot yield properly measured values.

NOTE

The automatic tracking function does not work at 3R and 38R, which are the boundaries of the corneal curvature radius measuring range. The measuring range boundary was measured manually using a model eye manufactured by HUVITZ.

NOTE

Measurement errors can occur under one of these conditions: out-of-focus, closed eyelid, tear film irregularity, high standard deviation in multiple measurements, movement, and measurement conducted outside of the range.


6.6. MAINTENANCE

6.6.1. After operation

1. Software exit and power off.


(1) On any screen, select the patient information icon () in the upper left corner.

(2) On the screen below, select the Power Off icon () in the bottom left corner.

(3) When you select the Power Off icon (), a pop-up window appears as shown below.


Selecting the Packing mode option will move the chin rest and the body to the lowest position before powering off.

NOTE

When placing the instrument in the packaging box, you must power it off by selecting the Packing mode and press OK.

2. If the monitor and/or other accessories are still connected, power them off.
3. Press the power off button (O) on the base plate to switch off the instrument.

6.6.2. Cleaning and Instrument Protection

There is a risk of electric shock if the HTG-1 device is not properly disconnected before cleaning the device. Therefore, ensure that the power plug is disconnected before cleaning the device.

1. Exterior

- (1) Ensure that the exterior is kept clean using a soft cloth. If the exterior has stains, wipe it with a soft cloth with neutral detergent diluted with water. When wiping the exterior, do not use volatile solutions, such as thinner or benzene.
- (2) Wipe the LCD monitor with a dry, soft cloth. Do not use a sponge or a wet cloth.
- (3) Do not press the LCD monitor with force or place a magnetic object near it.

2. Patient-contacted Areas

- (1) Wipe the headrest and chin rest clean where it makes a direct contact with patients, using clean gauze or cotton. If necessary, wet the cloth with rubbing alcohol and gently wipe the parts with it.
- (2) When the chin rest paper is in use, remove the used paper before examining the next patient.

3. Placido disk


- (1) The Placido disk dictates precision and is sensitive to pressure. The surfaces are susceptible to scratches.
- (2) Be very careful when cleaning the Placido's disk's surface. Use a dry, lint-free cloth when cleaning it.
- (3) If necessary, gently wipe the surface with a cloth slightly dampened with water.
- (4) Make sure that no dust gets into the small holes.

4. Others


- (1) If the instrument will not be used for an extended period, store it with a dust cover.
- (2) Before sending the instrument to an agency or HUVITZ for maintenance, clean the headrest and chin rest with rubbing alcohol.

6.7. Send Data to Web-Viewer


1. Fill in Server IP and Server Port information.
(The Server IP refers to the local IP address, where HIS S/W program is installed. The Server Port is fixed to 8080.)


2. Select the patients from the Patients List to send their data to the web viewer. Then, click the 'TRANSFER' button to send all of the selected patient data.

3. If only specific data need to be sent to the web viewer, select only those data from the data list.

4. When the data transfer is successful, the following message box will appear on the screen.

7

TROUBLESHOOTING GUIDE





If the device does not function properly, try the solutions listed in the table below before contacting the agency where you purchased the instrument.

If the device continues to malfunction after trying the solutions below, power it off and contact the agency.

The screen does not turn on.	The power cable may not be connected.	Check the connection status of the power cable.
	The power switch may not be turned on.	Check if the power switch is turned on.
The screen does not turn on, even though the power is on.	The instrument may be in Sleep mode.	Touching the screen will restore the system from Sleep mode.
The screen suddenly turns off.		
The image cannot be captured.	The patient may not be looking at the fixed target at the time of the measurement.	Instruct the patient to focus on the fixed target.
The measurement data is not visible.	The patient's eyelid or eyelashes may have interfered with the image clarity.	Instruct the patient to open their eyes wide. If the patient cannot open their eyes wide, lift the patient's eyelid while being careful not to apply pressure against their eyes.

Specifications and Accessories


8.1. Standard Accessories

 Chin rest paper	 Touch pen	 Power cable
 Spare Fuse (250V T 3.15A)	 Dust cover	 User manual
 Model eye	 Foot switch	

8.2. Specifications

Corneal Topography		
Working distance	80 ~ 100 mm	
Placido's disk	24 rings	
Keratometry		
Parameter	Measuring Range	Accuracy
Corneal measuring range	Up to 9.8 mm (R = 8.0 mm, 42.2 D, n = 1.3375)	-
Corneal curvature radius	3.00 to 38.00 mm (9.00 to 112.5 Diopter, n = 1.3375)	The accuracy and repeatability are in accordance with Type A, ISO 19980
Direction of principal meridians	0° to 180°	According to the ISO 10343:2014
White-to-white distance	7 to 14 mm	± 0.05 mm
Pupil diameter	0.5 to 10 mm	± 0.05 mm
Common		
Display	Swivel & Tilttable 10.1-inch, touch panel color LCD	
Horizontal movement	55 mm (back and forth), 100 mm (left and right)	
Vertical movement	30 mm	
Chin rest moving distance	62 mm (up and down), motorized	
Auto tracking	X, Y for positioning, Z for working distance	
Rated output	AC 100 - 240 V, 50/60 Hz, 1.6 - 0.7 A	
PC	Built-in computer	
Size	312(W) x 553(D) x 518(H) mm	
Weight	19 kg	

8.3. Drawings of System

9

EMC INFORMATION

Manufacturer announcement – electromagnetic waves trouble

• Electromagnetic waves trouble

HTG-1 should be used in the below mentioned electromagnetic wave environment. HTG-1 purchaser or user needs to confirm whether HTG-1 is used in this type of environment.

RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	Class A
Voltage fluctuations/flicker IEC 61000-3-3	Complies

• Electromagnetic waves tolerance

HTG-1 is to be used in the below designated electromagnetic wave environment. HTG-1 customer and user need to guarantee that the HTG-1 will be used in this type of environment.

Electrostatic discharge (ESD) IEC 61000 – 4 – 2	contact ± 8 kV in the air ± 15 kV	contact ± 8 kV in the air ± 15 kV
Electric rapid transients/burst IEC 61000 – 4 – 4	power supplying line ± 2 kV input/output line ± 1 kV	power supplying line ± 2 kV input/output line ± 1 kV
Surge IEC 61000 – 4 – 5	between lines ± 1 kV between line and grounding ± 2 kV	differential mode ± 1 kV common mode ± 2 kV
Voltage dip, instantaneous interruption, voltage fluctuation at the power input line IEC 61000 – 4 – 11	For 0.5 cycle < 5 %UT (UT's > 95 % decrease) For 5 cycle, 40 % UT (UT's 60 % decrease) For 25 cycle, 70 %UT (UT's 30 % decrease) For 5 seconds < 5 % UT (UT's > 95 % decrease)	For 0.5 cycle < 5 %UT (UT's > 95 % decrease) For 5 cycle, 40 % UT (UT's 60 % decrease) For 25 cycle, 70 %UT (UT's 30 % decrease) For 5 seconds, < 5 %UT (UT's > 95 % decrease)
Power frequency magnetic field (50/60 Hz) IEC 61000 – 4 – 8	30 A/m	30 A/m
Other UT is the A.C. power voltage for before approving the test level.		

- **Electromagnetic waves tolerance**

HTG-1 is to be used in the below mentioned electromagnetic wave environment. HTG-1 purchaser or user needs to confirm whether HTG-1 is used at this environment.

Conductivity RF electromagnetic field IEC 61000 – 4 – 6	3 Vrms 150 kHz~80 MHz	3 Vrms
Radioactivity RF electromagnetic field tolerance IEC 61000 – 4 – 3	10 V/m 80 MHz~2.7 GHz scope	10 V/m

SERVICE INFORMATION

Repair: If the problem has not resolved despite following the Troubleshooting guide in Chapter 7, please contact Huvitz's agent along with the following information:

- 1.1 Name of Equipment Type: Corneal Topographer HTG-1
- 1.2 Typical No. of Equipment: Typical number consisted of 8 digits and characters written on its name plate.
- 1.3 Explanation on its symptom: Description in detail.
- 1.4 Supply of parts required for repair:

* The preservation period of parts needed to repair this machine is seven (7) years after the product has been discontinued.

- 1.5 Parts to be repaired by qualified service manpower:

- 1.6 The following parts are consumable or may degrade in quality after extended use. Users should not attempt to replace them on their own. Please contact a Huvitz agent for replacement if these parts have been sufficiently consumed or have degraded due to prolonged use.

- 1.6.1 Back up battery for clock and data.

SOFTWARE LICENSE AGREEMENTS

The software including all the contents of files, disk or other media may be used only when you agree with the "SOFTWARE LICENSE AGREEMENT" below.

Before using the Software, thoroughly read "SOFTWARE LICENSE AGREEMENT"

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SOFTWARE LICENSE AGREEMENT

This SOFTWARE LICENSE AGREEMENT (hereafter referred to as Agreement) is a legal contract between you (regardless of private users or corporate users) and HUVITZ CO., LTD. (hereafter referred to as HUVITZ).

By installing, copying, downloading or in any way starting to use the Software, you accept and agree to all the terms and conditions of this Agreement. If you don't agree with the terms and conditions of this Agreement, return the disc package and the instruction manual to the location where you obtained them.

1. GRANT OF LICENSE

HUVITZ grants to you a non-exclusive, non-transferable limited license to use the Software under the terms and conditions of this Agreement.

- You may use the Software in a single location on a hard disk of a computer
- "Use" of the Software means reading the Software to a temporary memory such as RAM of a computer, or installing the Software on a fixed memory such as a hard disk.
- When the Software is used on multiple computers, you must obtain the same number of the licenses as the number of computers where the Software will be used even if the Software will not be used concurrently on multiple computers.
- Installing the Software into a fixed memory such as a hard disk of a computer from the server doesn't mean, "Using" the Software as long as the network server runs only for installing the Software into another computer. In such cases, once another computer has the software installed from the server, that computer needs the license.

2. INTELLECTUAL PROPERTY OWNERSHIP

You acknowledge that the Software and the Documentation and all the duplications are the intellectual property of HUVITZ and that HUVITZ owns any and all source codes, copyrights, trademarks, patents, trade secrets, confidential information and any other proprietary rights of the Software and the Documentation. The Software and the Documentation are protected by copyright laws and international copyright treaties, as well as other intellectual property laws and treaties.

3. SOFTWARE VERSION

The Software and the Documentation may be updated or modified at HUVITZ's own discretion from time to time without any prior notice.

4. RESTRICTIONS

- You may not use any HUVITZ trademarks, trade names or logos without the prior written consent of HUVITZ.
- You may not copy, duplicate, divert, modify, distribute, reverse engineer, decompose, disassemble or otherwise attempt to discover the source code of the Software and the Documentation without the prior written consent of HUVITZ.
- You may not use any HUVITZ trademarks, trade names or logos without the prior written consent of HUVITZ.
- You may not license or sublicense the Software, or assign, transfer or convey in any manner the Software to any third person without the prior written consent of HUVITZ. Any attempt to grant the Software to any third person will terminate this Agreement and relieve all obligations which may be undertaken on HUVITZ, regardless of whether the Software is paid for or not.

5. WARRANTY

- You may receive a replacement when you find any physical damage to the CD-ROM or printed materials within 30 days from the receipt of the Software.
- Unless otherwise provided herein, the Software is being delivered to you on an "AS IS" basis. ALL EXPRESS OR IMPLIED CONDITIONS, REPRESENTATIONS AND WARRANTIES, INCLUDING BUT NOT LIMITED TO ANY IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE OR NON-INFRINGEMENT OF THIRD PARTY RIGHTS ARE DISCLAIMED. HUVITZ WILL NOT WARRANT THAT THE OPERATION OF THE SOFTWARE IS UNINTERRUPTED, OR ERROR-FREE.

6. LIMITATION OF LIABILITY

IN NO EVENT WILL HUVITZ BE LIABLE TO YOU FOR ANY DAMAGES, LOSS, CLAIMS OR COSTS WHATSOEVER INCLUDING BUT NOT LIMITED TO ANY INCIDENTAL, INDIRECT, CONSEQUENTIAL OR SPECIAL DAMAGES, OR ANY LOST REVENUES, PROFITS, DATA OR INFORMATION, ARISING OUT OF OR RELATED TO THE USE OF OR INABILITY TO USE SOFTWARE, EVEN IF HUVITZ IS NOTIFIED OF THE POSSIBILITY OF SUCH DAMAGES, LOSS, CLAIMS OR COSTS.

7. FORCE MAJEURE

In no event HUVITZ will be liable to you for any delay or failure in the performance of its obligations under this Agreement if and to the extent such delay or failure in performance arise from any cause or causes beyond the reasonable control of HUVITZ, including, but not limited to, fire, explosion, flood, earthquake, war, rebellion or riots.

8. INDEMNIFICATION

You agree to defend, indemnify and hold harmless HUVITZ, its affiliates and their directors, employees and agents from any damages, loss, liabilities or expenses (including attorneys fee) resulting from, or incurred in connection with or based upon (a) your breach or non-compliance of this Agreement or (b) installation or other operations to set up the Software by persons other than HUVITZ personnel or those specified by HUVITZ.

9. TERM AND TERMINATION

This Agreement will commence when you open this package or start the use of the Software and continue in effect until terminated by either party. You may terminate this Agreement at any time for convenience and HUVITZ may terminate this Agreement at any time for convenience. This Agreement will be immediately terminated if you fail to comply with any term or condition in this Agreement. Your payment will not be refunded in any event of termination. When this Agreement is terminated, you will cease your use of Software, and destroy and erase all the Software and Documentation. Section 4, 5, 6, 7, 8, 13 and 14 will be unaffected and survive the termination of this Agreement.

10. NO ASSIGNMENT

This Agreement or any part of this Agreement may not be assigned or transferred by you without the prior written consent of HUVITZ. Any assignment or transfer without such consent will be null and void.

11. NO THIRD PARTY RIGHTS

This Agreement is intended to be solely for the benefit of the parties and is not intended to confer any benefits upon or create any rights in favor of, any person other than the parties in this Agreement.

12. EXPORT REGULATIONS

Software is subject to Japanese export security laws and regulations and may be subject to export laws and regulations in other countries. You agree to comply strictly with such export laws and regulations.

13. SEVERABILITY

If any Section in this Agreement is found to be void or unenforceable by applicable law, such Section will be deleted automatically and the remaining Sections of this Agreement will remain in full force and effect.

14. NO WAIVER

The delay, omission or failure of HUVITZ to enforce any right or remedy in this Agreement will not be construed to be a waiver of such right or remedy of HUVITZ.

15. ENTIRE AGREEMENT

This Agreement constitutes the entire agreement between the parties and supersedes any prior written or oral agreement between the parties concerning the Software.