

HD-18B

Features
Big table top ideal for placing larger equipment or two small equipment, like the auto refractometers.
Elegant and streamlined design;
Metal base makes it durable and long standing for use.
High quality casters provide safe and easy movement.



Specs
Table top size: 750mmx480mm; (29.5"x18.9")
Elevation range: 670mm-870mm; (26.4"x34.2")
Max Load: 50kg; (110LBS)
Voltage: 220V/50HZ or 110V/60HZ
Power: 50W
Weight: 20kg (44LBS)



HD-18C

Features
Two table top size available to meet the needs of installing different instrument at your disposal
Specially designed eyepiece holder built inside the drawer for convenient storage of additional eye piece and swop out in operation.



Specs
Table top size: 900x450mm;(35.4"x17.7")
1000x550mm;(39.4"x21.7")
Elevation range: 685mm-885mm; (26.9"-34.8")
Load: 50kg(110LB)
Voltage: 220V/50HZ or 110V/60HZ
Weight: 38kg (83.6LB) for small table top
42kg (92.4LB) for big table top

