


HD-40



Features
Long table top adopted to seat larger ophthalmic diagnostic equipment like the auto auto perimeters, OCT.
Wide range of elevation movement
Metal base makes it durable and long standing for use.
Heavy duty casters provide safe and easy movement.

Computer holding bracket underneath the table top is available upon the request.

Specs
Table top size: 920mmx575mm; (36.2"x22.6")
Elevation range: 665mm-865mm; (26.2"x34.1")
Max Load: 50kg; (110LBS)
Voltage: 220V/50HZ or 110V/60HZ
Power: 50W
Weight: 35kg (77LBS)



HD-40C



2 Instruments or 3 Instruments
Tabletop size: 1100x575mm
(43.3"x22.6")
Load capacity : 50kg (110lbs)
Elevation range: 660-860mm
(25.6"-33.8")
Optional computer cabinet and printer
holder available
Medical grade casters for faster and
quieter movement