

WCT-100

WCT-100 is a newly designed motorized instrument table specifically for the auto perimeter. More than that, specially designed metal base makes it wheel chair accessible, which facilitates a lot for the patient with disability to do the eye exam.
It comes with a printer holder as standard configuration.



WTC-50

WCT-50
1 Instrument
Tabletop size: 99x53.5mm (38.9"x21")
Load capacity : 50kg (110lbs)
Elevation range: 725-925mm (28.5"-36.4")
Wheelchair accessible
Medical-grade casters for quieter and faster movement



Specs
Elevation range: 700-900mm;
Tabletop size: 500x920mm
Printer holder: available as standard
Package dimensions: 1035x700x470mm
Net weight: 50kgs

