

Stool

Presenting Our Professional Series And Classic Series of Medical Stools

Hongdee proudly introduces our new Professional Series line of medical stools specifically designed for the Ophthalmology, Optometry, and dentistry markets. Our Professional Series combines with our Classic Series to offer a wide selection of medical stools for clinics, hospitals, and private practice eyecare and dental professionals. All of our stools are ergonomically designed with comfort, support, safety, and quality in mind. We are committed to the medical community we serve.



HS-100

With its 14.2” wide seat and 3.5” thick high-density Polyfoam cushion this stool offers a great deal of comfort for doctors and technicians. Easy height adjustment by the circular ring ranges from 20.8” to 28.1”, and 5-leg base stability,make HS-100 an ideal stool choice for healthcare professionals.

- Features
- Seating Comfort

14.2” diameter seat with 3.5” thick, high-density polyform provides excellent comfort

Up and down elevation Easy height adjustment using the continuous ring under the seat Effortlessly raise or lower the stool to the desired position.

5-Star aluminum Base:Robust and strong metal base with five legs to provide much more stability and even weight distribution.

Stability

25” aluminum base provides secure stable rolling and weight distribution.

Dimensions	(HS-100)
SEAT SIZE	36cm(14.2")
SEAT FOAM THICKNESS	9cm(3.5")
RETRACTED SEAT HEIGHT	53cm(20.8")
Extended seat height	71.5cm(28.1")
BACK SIZE	0
ELEVATION CONTROL	Ring lever
CASTER TYPE	PU
WEIGHT CAPACITY	150kgs(330lbs)
ELEVATION STROKE	18.5cm(7.3")
BACK ADJUSTMENT RANGE	/
Elevation range	53- 71.5cm(20.8"-28.1")
CYLINDER COLOR	Black
BASE COLOR	Black
CASTER COLOR	Black

