

HS-200 Pro

The HS-200 Pro Stool with Ring and Backrest is designed to provide optimum comfort and support during long-term seating when performing consecutive patient exams. Its 16” diameter round seat and 3.6” thick high-density polyfoam cushion promotes maximum sitting comfort.The adjustable backrest enhances spinal posture.Easily reach and grab anywhere on the ring without having to look to adjust the seat height up or down from a retracted height of 22.2” up to an extended height of 30”. This height adjustment allows doctors and technicians to keep their feet evenly grounded for better posture and blood flow.

Features

Sturdy 5-leg, 23.6” diameter polished aluminum alloy base provides greater stability and weight distribution up to 150kgs (330lbs) capacity.
2.5” polyurethane (PU) casters allow smooth rolling on carpet and various types of flooring.
Comfortable 16” wide, 3.6” thick poly-foam round seat for long-term cushioned sitting comfort.
Large 16.5” x 8.6” contoured backrest angles up and down, back and forth for comfortable spinal support, particularly during long sitting sessions.
360-deree circular ring within easy reach for up and down height adjustment from 22.2” to 30”.
Durable upholstery materials allow cleaning with alcohol cleansers for a safer hygienic environment without damaging the material.
CAL 133 Certified for fire code safety standards for upholstered furniture used in public buildings.

Model No	HS-200 Pro
Dimensions	Hongdee custom stool with back
SEAT SIZE	16""Wide (40cm)
SEAT FOAM THICKNESS	3.54" Thick (9cm)
SEAT SHAPE	Round
Foam Type	Poly Foam
BACK SIZE	42x22cm (16.5"x8.6"): Contoured
SEAT & Back Upholstery	Standard and alchol-proof Leather
SEAT & BACK COLORS	Black color
ELEVATION CONTROL	360-degree RING
CASTER TYPE	2.5" black caster
LEG TYPE	5-Star alloy base
BASE DIAMETER	23.6" (600mm) diameter
WEIGHT CAPACITY	330IBS
Retractable height	56.5cm(22.2")/18.7"(47.5cm)
Extended height	76.5cm(30")/24.2"(61.5cm)
ELEVATION STROKE	18.5cm (7.25")/5.3"(13.5cm)
Elevation range	56.5- 76.5cm(22.2" -30") /18.7" -24.2"(47.5-61.5cm)
WARRANTY	One year warranty
CYLINDER COLOR	Black
BASE COLOR	Black
CASTER COLOR	Black



HS-300

HS-300 comes with a 42x22cm (16.5”x8.5”) contoured back providing comfort at a great price. The modern look and ergonomically style seat (42.5x38.5cm (16.7”x15.2”) offer alot of value to eyecare and dental professionals. HS-300 is equipped with a tilt adjustable backrest that can be horizontally and vertically positioned for an ergonomic fit. The 5-leg robust aluminum and metal base gives it a stable and secure seating and moving experience.

Dimensions	(HS-300)
SEAT SIZE	42.5x38.5cm(16.7"x15")
SEAT FOAM THICKNESS	7cm(2.8")
RETRACTED SEAT HEIGHT	53cm(20.8")
Extended seat height	71.5cm(28.1")
BACK SIZE	42x22cm(16.5"x8.7")
ELEVATION CONTROL	Side lever
CASTER TYPE	PU
WEIGHT CAPACITY	150kgs(330lbs)
ELEVATION STROKE	18.5cm(7.3")
BACK ADJUSTMENT RANGE	13cm(5")
Elevation range	90-108cm(35.4" -42.5")
CYLINDER COLOR	Chrome
BASE COLOR	Chrome
CASTER COLOR	Black

