

HS-400

HS-400 is our stool model with a rounded backrest,height adjustment hand-grip lever, and a footrest ring for comfortable sitting with leg support. Several benefits of the stool footrest include reducing lower back and leg discomfort, allowing the feet to be elevated while sitting in an effortless neutral posture, promoting better blood circulation,preventing health concerns like varicose veins, and accommodating shorter individuals to sit more comfortably at a higher seat height while keeping their feet planted on the footrest.

Dimensions	(HS-400)
SEAT SIZE	36cm(14.2")
SEAT FOAM THICKNESS	7cm(2.8")
RETRACTED SEAT HEIGHT	53cm(20.8")
Extended seat height	71.5cm(28.1")
BACK SIZE	50x25cm(19.7"x9.8")
ELEVATION CONTROL	Side lever
CASTER TYPE	PU
WEIGHT CAPACITY	150kgs(330lbs)
ELEVATION STROKE	18.5cm(7.3")
BACK ADJUSTMENT RANGE	0
Elevation range	88-106cm(34.6"-41.7")
CYLINDER COLOR	Chrome
BASE COLOR	Chrome
CASTER COLOR	Black



HS-500

HS-500's 38x25cm (15" ×9.8") contoured backrest and 46.2x41.5cm (18.2" ×16.3") seat provides a pleasant seating experience. HS-500 is built with a tilt adjustable backrest that can be horizontally and vertically positioned for an ergonomic fit. The control lever under the seat allows easy up and down elevation. Additionally,HS-500 includes self-locking casters that automatically lock the stool into place when the user sits down, promoting stability and preventing unwanted chair movement during the exam.

Dimensions	(HS-500)
SEAT SIZE	46.2x41.5cm(18.2"x16.3")
SEAT FOAM THICKNESS	7cm(2.8")
RETRACTED SEAT HEIGHT	53cm(20.8")
Extended seat height	71.5cm(28.1")
BACK SIZE	38x25cm(15"x9.8")
ELEVATION CONTROL	Side lever
CASTER TYPE	PU
WEIGHT CAPACITY	150kgs(330lbs)
ELEVATION STROKE	18.5cm(7.3")
BACK ADJUSTMENT RANGE	13cm(5")
Elevation range	87-105cm (34.3"-41.3")
CYLINDER COLOR	Chrome
BASE COLOR	Chrome
CASTER COLOR	Black (Self-locking)

